



L2

COVID-19

Forensic Cleaning

FOR CONFIRMED COVID-19 CASES WITHIN YOUR WORKPLACE



HOTLINE: 1300 663 963



service@sps-cleaning.com



www.sps-cleaning.com



+ SERVICE SPECIFICATION

Level 2

In the event that a person is confirmed to have Covid-19 in your workplace, our team will perform a sterilising & biohazard deep disinfectant service for your asset. Our comprehensive program using STERI-7 XTRA* will effectively disinfect and eliminate 99.9999% of bacteria, viruses, and spores.

Sanitisation & disinfecting service specification:

- SPS cleaners wearing Biohazard Personal Protective Equipment (PPE)
- Broad-spectrum disinfection fogging service that targets all inaccessible areas
- Full disinfectant wipe of all frequently touched areas such as door handles, lift buttons, workstation desks, keyboards, mouse, phone, light switches, chairs and armrests
- Disinfectant mopping of hard floor areas
- Disinfectant wipe of all kitchen areas including microwaves, fridges etc.
- Discard all used PPE and any Biohazard waste on completion

FOR CONFIRMED CASES IN YOUR WORKPLACE

Seal off all areas where the confirmed case has visited

Step One: Risk Assessment

SPS staff will conduct a pre-start site risk assessment and job safety analysis worksheet. Our comprehensive risk assessment process identifies any potential hazards within the workplace and what protective control measures must be put in place before the commencement of any work. This ensures that all SPS employees and client personnel are as safe as possible at all times.

Step Two: Sterilising & Biohazard Clean

All SPS cleaners have gone through a specific Sterilising & Biohazard disinfectant cleaning training module. SPS will use our antimicrobial fogging service where available. However, this must take place after hours as Biohazard protective equipment must be worn by anyone present.

Step Three: Post-Clean Analysis

After the successful delivery of cleaning services, a post-clean analysis is conducted to ensure that the best result possible was achieved and to ensure the continuous improvement of service delivery. The post-clean analysis checks all cleaning methodologies used against the result achieved to ensure that all methods used for each task resulted in an excellent cleaning standard.

* If STERI-7 EXTRA is unavailable we will use a similar broad-spectrum disinfectant

STERI-7 XTRA Concentrate

Features and benefits

- Reactive Barrier Technology protection between cleans
- High-level disinfectant cleaner
- Non-corrosive
- Non-residual organoleptic effect on food
- Low toxicity
- Effective in soft or hard water
- No reported resistance
- Triple active reducing need to rotate products
- Maintains efficacy in heavy organic soiling, blood and proteins

Recommended usage

The STERI-7 XTRA Concentrate can be used wherever the highest standards of disinfection is required and are suitable for use in healthcare and food processing environments and all other workplaces where there is a risk of cross contamination. The product has been tested against and is effective against a number of commonly occurring bacteria, yeast and viruses that are known to be highly transmissible and can result in infections and illnesses.



HOTLINE
1300 663 963

EMAIL
service@sps-cleaning.com